

THUASNE®

# MOBIDERM

|    |   |    |
|----|---|----|
| fr | Plaque et bandage mobilisateurs pour les œdèmes .....                         | 3  |
| en | Mobilizing pad and bandage for edema .....                                    | 4  |
| de | Mobilisierende Pads und Bandagen zur Behandlung von Ödemen .....              | 6  |
| nl | Mobilisatiepad en -verband voor oedemen .....                                 | 7  |
| it | Placca e benda di mobilitizzazione per edemi .....                            | 9  |
| es | Placa y venda movilizadoras para edemas .....                                 | 10 |
| pt | Placa e ligadura de mobilização para os edemas .....                          | 12 |
| da | Mobiliserende plade og bind til ødemer .....                                  | 14 |
| fi | Mobilisoiva levy ja sidos turvotukseen .....                                  | 15 |
| sv | Mobiliserande ark och bandage för ödem .....                                  | 17 |
| el | Πλάκα και επίδεσμος κυνηγούντοινόντων για τα οδόντωτα .....                   | 18 |
| cs | Mobilizační podložka a bandáž pro léčbu edémů .....                           | 20 |
| pl | Arkusz i bandaż mobilizujący na obrzęki .....                                 | 21 |
| lv | Mobilizējošā plātnīte un pārsējs tūskām .....                                 | 23 |
| lt | Poodinė cirkuliacių skatinanti plokštelių ir juosta ištinusioms vietoms ..... | 24 |
| et | Nahaaluseid kudesid mobiliseeriv padi ja side ödeemide puhuks .....           | 26 |
| sl | Mobilizacijska blazinica in povoj za edeme .....                              | 27 |
| sk | Podložka a obväz na mobilizáciu edémov .....                                  | 29 |
| hu | Ödémacsökkentő mobilizáló lap és fáslí .....                                  | 30 |
| bg | Мобилизиращи подложки и бинтове за отоци .....                                | 32 |
| ro | Placă și bandaj de mobilizare pentru edeme .....                              | 34 |
| ru | Мобилизующие компрессы и повязки для лечения отеков .....                     | 35 |
| hr | Podložak i zavoj za imobilizaciju kod edema .....                             | 37 |
| zh | 用于水肿的活动泡沫片和活动绷带 .....   | 38 |
| ja | リンパ浮腫用活性化プレートおよび包帯 .....  | 39 |
| ko | 종 모빌라이징 패드 및 봉대 .....   | 41 |
| ar | لوح ولفافة تثبيت للوذمات .....  | 43 |

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en

## **MOBILIZING PAD AND BANDAGE FOR EDEMA**

#### **Description/Destination**

Mobiderm mobilizing devices are intended for the treatment of upper or lower limb edema.

Target population: adult patients and children suffering from indicated pathologies.

This range is available in several versions:

Mobiderm bandages:

- small blocks (5 x 5 mm) A
- big blocks (15 x 15 mm) B

Mobiderm pads:

- small blocks (5 x 5 mm) A
- big blocks (15 x 15 mm) B

The devices feature Mobiderm technology consisting of foam blocks encased between two sheets of non-woven material.

Mobiderm devices are available in different sizes of pads and bandages and in two foam block sizes to be used according to the needs and the body area treated.

#### **Composition**

Textile components: polyester - acrylic adhesive mass.

Foam components: polyurethane (large blocks) - polyethylene (small blocks) -

#### **Properties/Mode of action**

Mobiderm devices are intended to be used under a compression bandage. In this way, they provide subcutaneous drainage via Mobiderm technology, consisting of foam blocks positioned randomly between two sheets of non-woven material.

In contact with the skin, the blocks provide a pressure differential between the contact area beneath them and their periphery, creating multi-directional circulation channels for the lymph under the skin. The lymph travels from the areas pressed by the blocks towards the free areas. Efficient lymph drainage and circulation is thus promoted. C

#### **Indications**

Treatment of acute or chronic edema management (e.g: lymphedema).

#### **Contraindications**

Do not use in case of major dermatological conditions on an affected limb.

Do not use in the event of known allergy to any of the components.

Do not use in case of severe peripheral neuropathy of the affected limb.

Do not use in the event of septic thrombosis.

Do not use in the event of phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of an affected limb.

Specific to the lower limbs:

Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) < 0.6.

Do not use in the event of decompensated heart failure.

Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg). Do not use in case of extra-anatomical bypass of the affected limb.

Specific to the upper limbs:

Do not use in case of brachial plexus condition. Do not use in case of vasculitis of the extremities.

#### **Precautions**

Strictly comply with your healthcare professional's prescription and recommendations for use.

The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

In the event of discomfort, significant restriction, pain, variation in limb volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product, remove the product and seek the advice of a healthcare professional.

Do not apply the product directly to an open wound without a dressing.

Application of this device can cause the appearance of marks on the skin. These marks are benign and temporary; do not stop treatment without medical advice.

Due to the constant evolution of the pathology, regular follow-up with a healthcare professional is recommended to check that the compression device is still suitable for the pathology.

Do not use when swimming.

Do not use products on the skin (creams, ointments, oils, gels, patches...) before applying the device, it could damage the product.

In the event of associated advanced venous disease or venous ulcer, seek the advice of a healthcare professional for appropriate compression treatment on the entire leg.

For hygiene, security and performance reasons, do not re-use the product for another patient.

#### **Undesirable side-effect**

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness. Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

#### **Instructions for use/Application**

It is recommended that your device be washed before it is used for the first time.

Verify the product's integrity before every use. A cotton bandage or protective bandage must be placed underneath the Mobiderm device.

It is possible to place cushioning devices (such as Varico, foam bandage) underneath Mobiderm devices.

It is possible to cut the devices to fit the area to be treated according to needs.

It is recommended that the devices be cut with scissors between 2 rows of foam blocks.

Check that there are no folds in the device once in place.

#### **Application of Mobiderm bandages:**

Apply the bandage circularly along the length of the limb without pressure and with an overlap of one row of blocks between each turn.

On the joints, we recommend an overlap of 2 rows of foam blocks in order to maintain the position of the device during various body movements.

Cut off surplus bandage if necessary.

Then apply the compression bandage according to the manufacturer's recommendations.

#### **Application of Mobiderm pads:**

Locate the area to be treated then cut the pad according to the shape of the chosen surface. Apply the pad to the area to be treated without pressure.

Then apply the compression bandage according to the manufacturer's recommendations.

#### **Care/Maintenance**

To ensure your Mobiderm devices last as long as possible and to prevent infection, it is important to wash them every day, or after each use.

Soak the device in tepid soapy water, rub with a clean sponge and squeeze out excess water.

Do not use chlorine-based detergents or harsh products. Air dry, away from direct heat (e.g. radiator, sun, etc.). Do not tumble-dry. Do not iron.

#### **Storage**

Store at room temperature, preferably in the original packaging.

#### **Disposal**

Dispose in accordance with local regulations.

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부위에서 공간이 있는 부위로 이동합니다. 그 결과 배수 효율과 림프 순환 측정됩니다. ©
THUASNE UK Ltd**
Thuasne SAS – SIREN/RCS Nanterre 542 091 186